

For Research Use Only

TGF beta Receptor I Polyclonal ANTIBODY



Catalog Number: 55391-1-AP

3 Publications

Basic Information

Catalog Number:
55391-1-AP

Size:
48 µg/150 µl

Source:
Rabbit

Isotype:
IgG

Purification Method:
Antigen affinity purification

Immunogen Catalog Number:

GenBank Accession Number:
NM_004612

GeneID (NCBI):
7046

Full Name:
transforming growth factor, beta receptor 1

Calculated MW:
56 kDa

Observed MW:
56 kDa

Recommended Dilutions:

Applications

Tested Applications:
ELISA

Cited Applications:
IHC, WB

Species Specificity:
human, mouse

Cited Species:
human, mouse

Note: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Na Tang	28942050	Exp Parasitol	WB
Jingjing Wu	29077259	J Cell Mol Med	WB
Boshi Wang	30745168	EBioMedicine	WB, IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data