

For Research Use Only

GSTO1 Monoclonal antibody

Catalog Number: 60150-1-Ig



Basic Information

Catalog Number:	60150-1-Ig	GenBank Accession Number:	BC000127	Purification Method:	Protein A purification
Size:	150ul	GeneID (NCBI):	9446	CloneNo.:	1E3E9
Source:	Mouse	UNIPROT ID:	P78417		
Isotype:	IgM	Full Name:	glutathione S-transferase omega 1		
Immunogen Catalog Number:	AG7199	Calculated MW:	28 kDa		

Applications

Tested Applications:
ELISA
Species Specificity:
human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

