

For Research Use Only

Anti-Mouse CD24 (M1/69)

Catalog Number: 65121-1-Ig



Basic Information

Catalog Number:

65121-1-Ig

Size:

100ug , 500 µg/ml

Source:

Rat

Isotype:

IgG2b, kappa

GenBank Accession Number:

BC075622

GeneID (NCBI):

12484

Full Name:

CD24a antigen

Purification Method:

Affinity purification

CloneNo.:

M1/69

Recommended Dilutions:

FC: 0.5 ug per 10⁶ cells in 100 µl suspension

Applications

Tested Applications:

FC

Species Specificity:

mouse

Positive Controls:

FC : mouse splenocytes,

Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ mouse splenocytes were surface stained with 0.5 ug Anti-Mouse CD24 (M1/69) (65121-1-Ig, Clone: M1/69) or 0.5 ug Rat IgG2b Isotype Control (LTF-2) (65211-1-Ig, Clone: LTF-2), and FITC anti-rat IgG2b Antibody. Cells were not fixed.