

For Research Use Only

HRAS-Specific Monoclonal ANTIBODY



Catalog Number: 66199-1-Ig

Basic Information

Catalog Number:

66199-1-Ig

Size:

150UL, Concentration: 740 µg/ml by Bradford method using BSA as the standard;

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC006499

GeneID (NCBI):

3265

Full Name:

v-Ha-ras Harvey rat sarcoma viral oncogene homolog

Calculated MW:

21 kDa

Purification Method:

Protein G purification

CloneNo.:

1B7E4

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data