

For Research Use Only

APC Anti-Human CD1c (L161)

Catalog Number: APC-65176



Basic Information

Catalog Number:

APC-65176

Size:

100 tests, 10 µl/test

Source:

Mouse

Isotype:

Mouse IgG1 kappa

GenBank Accession Number:

BC126465

GeneID (NCBI):

911

Full Name:

CD1c molecule

Purification Method:

The purified antibody is conjugated with allophycocyanin (APC) under optimum conditions. The conjugate is purified by size-exclusion chromatography.

CloneNo.:

L161

Applications

Tested Applications:

FC

Species Specificity:

Human

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light.

Storage Buffer:

PBS with 0.1% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human peripheral blood lymphocytes were surface stained with PE Anti-Human CD19 (PE-65110, Clone:H1B19) and 10.00 ul/test APC Anti-Human CD1c (APC-65176, Clone:L161). Cells were not fixed.



1X10⁶ human peripheral blood lymphocytes were surface stained with PE Anti-Human CD19 (PE-65110, Clone:H1B19) and APC-mouse IgG1 isotype control antibody. Cells were not fixed.