

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.



RECOMBINANT HUMAN TMPRSS3

Basic Information

Catalog Number:
Ag0443

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
E. coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
MGENDPPAVEAPFSFRSLFGLDDLKISPVAPDADAVAAQILSLLP
LKFFPIIVIGIILALALAIQLGIHFDCSGKYRCRSSFKCIELIARCD
GVSDCKDGEDEYRCVIRVGGQNAVLMQVFTAASWKTMCSDDWK
GHYANVACAQLGFPSYVSSDNLRVSSLEGQFREFVSDHLLP
DDKVTALHHSVYVREGCASGHVTLQCTACGHRRGYSSRIVG
GNMBSLSSQWVWQASLQFGYHLCGGSVITPLWIIAAHCYVDL
YLPKSWTIQVGLVSLDNPAPSHLVEKIVYHSKYKPKRLGNDIAL
MKLAGPLTFNEMQPVCLPNSEENFPDGVKCVTSGWGATEDG
GDASPVLNHAAVPLISNKICNHRDYYGGIISPSMLCAGYLTGGV
DSCQGDGGPLVCQERRLWKLVGATSGFGCAEVNKPGVYTR
VTSFLDWIHEQMERDLKT
(N-term-500aa encoded by BC074846)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),
or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.