

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RECOMBINANT HUMAN XIAP

Basic Information

Catalog Number:
Ag0725

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
E. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
MTFNSEFGSKTCVPADINKEEEFVEEFNRLKTFANFPSPGSPVSA
STLARAGFLYTGEGDTRCFSCHAAVDRWQYGDSAVGRHRK/
SPNCRFINGFYLENSATQSTNSGIQNGQYKVENYLGSRDHFAL
DRPSETHADYLLRTGQVVDISDTIYPRNPAMYSEEARLKSFQV
VFPDYAHLTPRELASAGLYYTGIDQVQCFCCGGKLNWEPED
RAWSEHRRHFPCFFVLRNLNIRSEDAVSSDRNFPNSTNL
PRNPSMADYEARIFFTGTWYISVNIKEQLARAGFYALGEGDKVK
CFHCGGGLTDWKPSEDPWEQHAHWYPGCKYLLEQKQOEYIN
NIHLTHSLEECLVRTTEKTPSLTRRIDDTIFQNPMMGEAIRMGFSF
KDIKIMEEKIQSGSNYKSLELVADLVNAQKDSMQDESSQTSL
QKEISTEEQLRRLQEEKLCIKMDRNIAIVFPCGHLVTCQCAE
AVDKCPMCTVITFKQIFMS
(1-497 aa encoded by BC032729)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),
or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.