

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RECOMBINANT HUMAN AGPAT6

Basic Information

Catalog Number:
Ag10294

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
E. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
KEFM~~SK~~HVHLMCYRICVRLTAITYHDREN~~R~~FRNGGICVANHT
SPIDVILASDGY~~Y~~AMVGVHGGLMGMQRAMKACPHWFER
SEVKDRHLVAKRLTEHVQDKSKLPI~~L~~IFPEGTCINNTSMVFKKG
SFEIGATVY~~P~~VAIKYDPQFGDAFWNSSKYGM/TYLLRMMITSWA
IVCSWY~~L~~PPMITREADEDAVQFANRVKSAIARQGGGLVDLLWDG
GLKREKVKDTFKEEQQLYSKIMVGNHKDRSRS
(206-456 aa encoded by BC051377)

Selected Validation Data

For technical support and original validation data for this product, please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),
or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.