

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RECOMBINANT HUMAN FIGNL1



Basic Information

Catalog Number:
Ag11709

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
E. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
KSLGASRSRGI L G F V P P I P K Q D G G E Q N G G M Q C K P Y G A G P T E P
A H P V D E R L K N L E P K M E L I M N E I M D H G P P M N W E D I A G V E F A K A T I
K E I V W P M L R P D I F T G L R G P P K G I L L F G P P G T G K T L I G K C I A S Q S G
A T F F S I S A S S L T S K W V G E G E K M R A L F A V A R C Q Q P A M I F I D E I D S
L L S Q R G D G E H E S S R R I K T E F L V Q L D G A T T S S E D R I L V G A T N R P
Q E I D E A A R R R L V K R L Y I P L P E A S A R K Q I V I N L M S K E Q C C L S E E E E
Q I V Q Q S D A F S G A D M T Q L C R E A S L G P I R S L Q T A D I A T I T P D Q V R P I
A Y I D F E N A F R T V R P S V S P K D L E L Y E N W N K T F G C G K
(326-674 aa encoded by BC051867)

Selected Validation Data

For technical support and original validation data for this product, please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),
or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.