

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human BHLHE40 protein



Basic Information

Catalog Number:
Ag12236

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, T-HIS, with N-terminal 6*His.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Peptide Sequence:
MERIPSAQPPPACLPKAPGLEHGDLPGMYP AHMYQV
YKSRRGIKRSEDSKETYLPHRLIEKRRDRINECIAQLK
DLLPEHLKLTTLGHLEKAVVLELTKHV KALTNLIDQQQ
QKIALQSGLQAGELSGRNVETGQEMFCSGFQTCAR
EVLQYLAKHENTRDLKSSQLVTHLHRVVSELLQGGTS
RKPSDPAPKVMDFEKEKPSPAKGESEGPKNCPVPIQR
TFAHSSGEQSGSDTDTDSGYGGESEKGDLRSEQPCF
KSDHGRRFTMGERIGAIKQESEEPPTKKNRMQLSDDE
GHFTSSDLISSPFLGPHPHQPPFCLPFYLIPPSATAYLP
MLEKCYWYPTSPVLYPGLNASAAALSSFMNPKISAP
LLMPQRLPSPLPAHPSVDSSVLLQALKPIPLNLETKD
(1-412 aa encoded by BC082238)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W: ptglab.com
(toll free in USA), or 1(312) 455-8498
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.