

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RECOMBINANT HUMAN ACP1



Basic Information

Catalog Number:
Ag17539

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
E. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:

GNI~~C~~RSP~~I~~AEAVFRKLVTDQNI~~S~~ENW~~M~~DSGAVSDW~~N~~VGRSPD
PRAV~~S~~CLRNIHG~~I~~HTA~~H~~KARQITKEDFATFDYILCMDES~~N~~L~~R~~DLN
R
(15-102 aa encoded by BC007422)

Selected Validation Data

For technical support and original validation data for this product, please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),
or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.