

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RECOMBINANT HUMAN JAK1



Basic Information

Catalog Number:

Ag18180

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

E. coli-derived, PET28a, with N-terminal 6^{His}.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

EDCNAMAFCAKMRSSKITEVNLEAPEPGVEVIFYLSD
REPLRLGSGEYTAEEELCIRAAQACRISPLCHNLFALYDE
NTKLWYAPNRTITVDDKMSLRHLHYRMRFYFTNWHGT
NDNEQSVWRHSPKKQKNGYEKKKIPDATPLLDASSLE
YLFAQQQYDLVKCLAPIRDPKTEQDGHDIENECLGMA
VLAISHYAMMKMQLPPEPKDISYKRYIPETLNKSIQR
NLLTRMRINNVFKDFLKEFNKTI CDSSVSTHDLKVKY
LATLETLTKHYGAEIFETSMLLISSENMNWFHSNDGG
NVLYYEVMTGNLGIQWRHKPNVVSVEKEKNLKRKK
(8-363 aa encoded by BC132729)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),
or 1(312) 455-8498 (outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.