

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RECOMBINANT HUMAN CDC14B



Basic Information

Catalog Number:

Ag18227

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

E. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

PGSVIGPQQFLVMKQTNLWLEGDYFRQKLGQEN
GQHRAAFSKLLSGVDDISINGVENQDQEQEPEPYSDD
DEINGVTQGDRLRALKSRRQSKTNAIPLTVILQSSVQS
CKTSEPNI S G S A G I T K R T T R S A S R I K S S V K S L S I S R T I K T V L
R

(349-498 aa encoded by BC156666)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),
or 1(312) 455-8498 (outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.