

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human EML4 protein



Basic Information

Catalog Number:

Ag18463

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

GLSDDSISAASTSDVQDRLSALESRVQQQEDEITVLK
AALADVLRRLAISEDHVASVKKSVSCKGQPSPRVIVP
MSCITNGSGANRKPSTSAVSIAGKETLSSAAKSGTE
KKKEKPGQGQREKKEESHNDQSPQIRASPSQPSSQP
LQIHRQTPESKNATPTKSIKRPSPAEKSHNSWENSDD
SRNLSKIPSTPKLIPKVTIKTADKHKDVIINQEGEYIKM
FMRGRPITMFI PSDVDNYDDIRTELPPEKLEWAYGY
RGKDCRANVYLLPTGEIVYFIASVVVLFNYEERTQRHY
LGHTDCVKCLAIHPDKIRIATGQIAGVDKDRPLQPH
VRVWDSVTLTLQIIGLGTFERGVGCLDFSKADSGVH
LCVIDDSNEHMLTVWDWQKKAKGAEIKTTNEVVLAV
EFHPTDANTIITCGKSHIFFWTWSGNLTKRGIFGKY
EKPKFVQCLAFNGDVLTDGDSGGVMLIWSKTTVEPT
PGKGPKVYQISKQIKAHDGVSFTLCQMRNGMLLTG
GGKDRKIIWLDHDLNPEREIEVPDQYGTIRAVAEGKA
DQFLVGTSRNFILRGTNDGFQIEVQGHDELWGLAT
HPFKDLLTCAQDRQVCLWNSMEHRLEWTRLVDEPG
HCADFHPSGTVVAIGTHSGRWFVLDAAE
(6-663 aa encoded by BC140845)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.