

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.



# RECOMBINANT HUMAN PAN-PAX

## Basic Information

**Catalog Number:**  
Ag19007

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*E. coli*-derived, PGEX4T, with N-terminal GST.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Validated Application:**  
Blocking peptide

**Peptide Sequence:**  
MPHNSIRSGHGGLNQLGGAFVNGRPLPEVWRQRIMDLAHQGV  
RPCDISRQLRVSHGCVSKILGRYYETGSIIRPGVIGGSKPKVATP  
KWVEKIGDYKRQNPIMFAWEIRDRLLAEGVCDNDTVPSVSSIN  
RIIRTKVQQPFNLPMDSCVATKSLSPGHTLIPSSAVTPPESPQSD  
SLGSTYSINGLLGIAQPGSDKRMDDSDQDSCRLLSIDSGSSSSG  
PRKHRLRTDAFSQHILEPLECPFERQHYPEAYASPSHTKGEQGL  
YPLPLLNSTLDDGKATLTPSNTPLGRNLSTHQTYPWADPHSPF  
AIKQETPEVSSSSSTPSSLSAFLDLQQVGSVPPFNAPPHAAS  
VYGGFTGQALLSGREIM/GPTLPGYPPHIPTSGQGSYASSAIGM  
VAGSEYSGNAYGHIPYSSYSEAWRFPNSSLSSPYYSSTSR  
PSAFPPTTATAFDHL  
(1-450 aa encoded by BC001060)

## Selected Validation Data

For technical support and original validation data for this product, please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),  
or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.