

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RECOMBINANT HUMAN JNK



Basic Information

Catalog Number:
Ag21426

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
E. coli-derived, PET30a, with N-terminal 6*His.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
MRSKRDNINFYSVEIGDSTFTLTKRYQNLKPIGSGAQQIVCAAY
DAILERNVAIKLSRPFQNTTHAKRAYRELVLMKCVNHKNIIQLL
NVFTPQKSLLEEFQDVYIVMELMDANLCQVQVMDLHERMSYLLY
QMLCGIKHLHSAGIHRDLKPSNIWKSDDLKILDFGLARTAGTS
FMIMTPYVTRYYRAPEVLGMGYKENVDIWSVGCIMGEMKGG
VLFPGTDHIDQWNVIEQLGTPCPEFMKLLQPTVRYVENRPKY
AGYSFEKLFDPVLPADSEHNKLSAQARDLLSKMLVIDASKRIS
VDEALQHPYINWVYDPSEAEAPPKIPDKQLDEREHTIEEWKEL
IYKEVMDLEERTKNGVIRGQPSPLGAAMNGSQHPSSSSVNDV
SSMSTDPFLASDTSLEAAAGFLGCCR
(1-427 aa encoded by BC130572)

Reconstitution and Storage

Reconstitution:
Reconstitute at 0.25 µg/µl in 200 µl sterile water for short-term storage.
Reconstitution with 200 µl 50% glycerol solution is recommended for longer term storage (see Stability and Storage for more details).
If a different concentration is needed for your purposes please adjust the reconstitution volume as required (please note: the ion concentration of the final solution will vary according to the volume used).
Note: Centrifuge vial before opening. When reconstituting, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution.

Shipping:
The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature (see below).

Stability and Storage


Store for up to 12 months at -20°C to -80°C as lyophilized powder.

Storage of Reconstituted Protein

Short Term Storage:
Store at 2-8°C for (1-2 weeks).

Long Term Storage:
Aliquot and store at -20°C to -80°C for up to 3 months, buffer containing 50% glycerol is recommended for reconstitution. Avoid repeat freeze-thaw cycles.

Selected Validation Data



For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),
or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.