

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.



# RECOMBINANT HUMAN PDE7A

## Basic Information

Catalog Number:  
Ag24334

Size:  
50 µg

Form:  
Available lyophilized

Species:  
human

Expression Source:  
*E. coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:  
Not tested

Endotoxin Level:  
Please contact the lab for more information

Validated Application:  
Blocking peptide

Peptide Sequence:  
MGI~~T~~L~~I~~W~~C~~L~~A~~L~~V~~L~~I~~K~~W~~I~~T~~S~~K~~R~~R~~G~~A~~I~~S~~Y~~D~~S~~S~~D~~Q~~T~~A~~L~~Y~~I~~R~~M~~L~~G~~D~~V~~R~~V  
R~~S~~R~~A~~G~~F~~E~~S~~E~~R~~  
(1-55 aa encoded by BC126360)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),  
or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.