

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.



RECOMBINANT HUMAN SCRIB

Basic Information

Catalog Number:
Ag25795

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
E. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
HPPEDAPAQP PTPGPAASPE QLSFRERQKY FELEVRVPGA E
GPPKRVSLV GADDLRKMQE EEARLQCKR AQMLREAAEA GA
EARLALDG ETLGEEEQED EQPPWASPPS TSRQSPASPP PLG
GGAPVRT AKAERRHQER LRVQSPEPPA PERALSPAEL RALE
AEKRAL WRAARMKSLQ QDALRAQML SRSQEGRGTR GPLL
RLAEAP SPAPTSPPTP VEDLGPQTST SPGRLSPDFA EELRSLE
PSP SPGQEEEDGE VALVLLGRPS PGAVGPEDVA LCSSRRPV
RP GRRGLGPVPS
(1331-1630 aa encoded by BC036905)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),
or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.