

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# RECOMBINANT HUMAN TRIM45



## Basic Information

**Catalog Number:**

Ag28300

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*E. coli*-derived, PET28a, with N-terminal 6<sup>+</sup>His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

VCQDCVVGEHREHPCDFTSNVIHKHGDSVWELLKGT  
QPHVEALEEALAQIHIIINSALQKRVEAVAADVRTFSE  
GYIKAIIEHRDKLLKQLEDIRAQKENSLLQKQLEQL  
LADMRTGVFEHLLTSGSDLEILITKRVVVERLRKLNK  
VQYSTRPGVNDKIRFCPQEKAGQCRGYEIDLHRARE  
KQTASFLLCKDAAGEIMGRGGDNVQVAVVPKDKKD  
SPVRTMVQDNKDGTYIISYTPKEPGVYTVVWCIKEQ  
HVQGSPTVTVRRKHRPHSGVFNCTFCSSGGQKTA  
RCACGGTMPGGYLGCGHGHKHPGHPHWSCCGKF  
NEKSECTWTGGQSAPRSLLRVAL  
(210-562 aa encoded by BC034943)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com) W:  
[ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.