

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RECOMBINANT HUMAN ANGPTL4



Basic Information

Catalog Number:

Ag30070

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

E. coli-derived, PET28a, with N-terminal 6⁺His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MSGAPTAGAALMLCAATAVLLSAQGGPVQSKSPRFA
SWDEMNVLAHGLLQLGQGLREHAERTRSQLSALERRL
SACGSACQGTGEGSTDLPLAPESRVDPEVLHSLQTQLK
AQNSRIQQLFHKVAQQQRHLEKQHLRIQHLSQSFGL
LDHKHLDHEVAKPARRKRLPEMAQPVDPAHNVSRLH
RLPRDCQELFQVGERQSGLFQIQPGSPPFLVNCMT
SDGGWTVIQRRHDGSDVDFNRPWEAYKAGFGDPHGE
FWLGLEKVHSITGDRNSRLAVQLRDWDGNAELLQFS
VHLGGEDTAYSLLTAPVAGQLGATTVPSSGLSVPFS
TWDQDHDLRDKNCAKSLSGGWWFGTCSHSNLNG
QYFRSIPQQRQKL

(1-374 aa encoded by BC023647)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.