

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.



RECOMBINANT HUMAN AGFG1

Basic Information

Catalog Number:
Ag3352

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
E. coli-derived, PKG, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
SAKRKQEEKHLKMLRDMTGLPHNRKCFDCDQRGPTYVNMIVG
SFVCTSCSGSLRGLNPPHRVKSISMTTFQQEIEFLQKHGNEVC
KQIWLGLFDDRSSAIPDFRDPQVKKEFLQEKYEKRWYVPPEQ
AKVWASVHASISGSSASSTSTPEVKPLKSLLGDSAPTLHLNKGT
PSQSPVWGRSQGQQEKQFDLLSDLGSDIFAAPAPQSTATAN
FANFAHFNSHAAQNSANADFANFDAGQSSGSSNFGGFPTAS
HSPFPQPQTGGSAASVNAVFAHFDNFFKSSADFQTFNTSQ
(4-303 aa encoded by BC030592)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),
or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.