

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# RECOMBINANT HUMAN ALS2



## Basic Information

**Catalog Number:**  
Ag5087

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*E. coli*-derived, T-HIS, with N-terminal 6\*His.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Validated Application:**  
Blocking peptide

**Peptide Sequence:**  
MDSKIKRSSTAEAGSKERGLVHIWQAGSFPITPERLPQWGGKTV  
LQAALGVKHGVLITDGEVYSFGTLPWRSQVPCSSPILENAL  
VGQYVITVATGSEFHSQAVTDNGVAYMCGENSAGQCAVANQQY  
VPEFNPVSIADSEASPLAVRILQLACGEEHTLALSISREIWAQW  
TGCCQLGLITTAFFVTKPQKVEHLAGRVLQVACGAFHSLALVQCL  
PSQDLKPVPERCNQCSQLLITMDKEDHVISDSHCCPLGVLTLE  
SQAENHASTALSPSTETLDRQEEVFENTLVANDQSVATE  
(1-304 aa encoded by BC029174)

## Selected Validation Data

For technical support and original validation data for this product, please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.