

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# RECOMBINANT HUMAN ADIPOR2

## Basic Information

Catalog Number:  
Ag5744

Size:  
50 µg

Form:  
Available lyophilized

Species:  
human

Expression Source:  
*E. coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:  
Not tested

Endotoxin Level:  
Please contact the lab for more information

Validated Application:  
Blocking peptide

Peptide Sequence:

MNEPTENRLGCSRTEPEPDIRLRKKGHLQDGTTRRGDNDSHQGDLEPILEASVLSHHKKSSEEHEY SDEAPQEDGFMGMSPLLQAH  
HAMEKVEEFVCKWEGRWVRMPHDVLPDWLKDNDFLHGH  
PPMPSFRACFKSIFRIHTETG  
(1-147 aa encoded by BC051858)

## Selected Validation Data

For technical support and original validation data for this product, please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),  
or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.