

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RECOMBINANT HUMAN AFF4



Basic Information

Catalog Number:

Ag6336

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

E. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MNREDRNVLRMKERERRNQEIQQGEDAFPPSSPLFAE
PYKVTSKEDKLSRIQSMGLGNYDEMKDFIGDRSIPKLV
AIPKPTVPPSADEKSNPNFFEQRHGGSHQSSKWTPV
GPAPSTSQSQKRSSGLQSGHSSQRTSAGSSSGTNSS
GQRHDRESYNNSSGSSSRKKGQHGSEHSKSRSSSPGKP
QAVSSLNSSHSRSHGNDHHSKEHQRSKSPRDPDAN
WDSPSRVPFSSGQHSTQSFPPSLMSKSNMLQKPTA
YVRPMDGQESMEPKLSEHYSSQSHGNSMTELKPSK
AHLTKLQIPSQPLDASASGDVSCVDEILKEMTHSWPP
PLTAIHTPC

(1-338 aa encoded by BC063007)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),
or 1(312) 455-8498 (outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.