

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RECOMBINANT HUMAN CONNEXIN 43



Basic Information

Catalog Number:
Ag6698

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
E. coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
MGDWSALGKLLDKVQAYSTAGGKWWLSVLFIFRILLGLTAVESA
VGD EQSAFR CNTQQPGCENV CYDKSFPI SHVRFWMLQIFVSV
PTLLYL AHV F YVMRKEEKL NKKEELKVAQTDGVMNDMHLKQIEI
KFKY GIEEHGKVMRGGLLR TYIIISILFKSIFEVAFLLIQWYIYG
FSLSAVYTCKRDPCHQMD CFLSRPTEKTIIFIMLVSLVSLALN
IIELFYVFFKGVKDRVKGKSDPYHATSGALSPAKDCGSQKYAYF
NGCSSPTAFLSPVSPPGYKLVTDGRNNSSCRNYNKQASEQNW
ANYSAEQNRMGQAGSTISNSHAQPFDFPDDNQNSKLAAGHE
LQPLAIVDQRFPSSRASSRASSRPRPDDLEI
(1-382 aa encoded by BC026329)

Selected Validation Data

For technical support and original validation data for this product, please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),
or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.