

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RECOMBINANT HUMAN PRKAB1



Basic Information

Catalog Number:
Ag7707

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
E. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
MGN TSSERAALERHGGHKTPRRDSSGGTKDGD RPKILMDSPE
DADLFHSEEIKAPEKEEF L AWQHDL E V NDKAP A Q A R P T V F R W T
GGGKEVYLSGSFNWWSKLP LTRSHN FVA I L D L P E G E H Q Y K F F
VDGQWTHDPSEPIVTSQLGT V N N I I Q V K K T D F E V F D A L M D S Q K
CSDVSELSSSPGPYHQEPYVCKPEERFRAPPILP P H L L Q M L N
KDTG I S C D P A L L P E P N H V M L N H L Y A L S I K D G V M M L S A T H R Y K K K
YV T T L L Y K P I
(1-270 aa encoded by BC001056)

Selected Validation Data

For technical support and original validation data for this product, please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),
or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.