

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# RECOMBINANT HUMAN DUSP22



## Basic Information

Catalog Number:  
Ag9831

Size:  
50 µg

Form:  
Available lyophilized

Species:  
human

Expression Source:  
*E. coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:  
Not tested

Endotoxin Level:  
Please contact the lab for more information

Validated Application:  
Blocking peptide

Peptide Sequence:

MGNGMNKILPGLYIGNFKDARDAEQLSKNKVTHILSVHDSARFV  
LEGVKYLCIPAADSPSQNLTRHFKEIKFHECRLRGESCLVHCL  
AGVSRSVTLVIAYIMTVIDFGWEDALHTVRAGRSCANPNVGFQ  
RQLQEFEKHEVHQYRQVLMKEEYGESPLQDAEEAKNILAAPGL  
KFWAFLRRL

(1-184 aa encoded by BC022847)

## Selected Validation Data

For technical support and original validation data for this product, please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),  
or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.