

For Research Use Only

CoraLite®488-conjugated TDP-43 (for IF/FC) Recombinant antibody



Catalog Number: **CL488-80002**

Basic Information

Catalog Number: CL488-80002	GenBank Accession Number: BC001487	Purification Method: Protein A purification
Size: 100ul , Concentration: 500 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 23435	CloneNo.: 16A22
Source: Rabbit	Full Name: TAR DNA binding protein	Recommended Dilutions: IF 1:50-1:500
Isotype: IgG	Calculated MW: 43 kDa	Excitation/Emission maxima wavelengths: 488 nm/515 nm
Immunogen Catalog Number: AG1231	Observed MW: 43 kDa	

Applications

Tested Applications: FC, IF	Positive Controls: IF : mouse brain tissue,
Species Specificity: Human, mouse, rat	

Background Information

The TARDBP gene encodes the TDP-43 protein, initially found to repress HIV-1 transcription by binding TAR DNA. TDP-43 has since been shown to bind RNA as well as DNA, and have multiple functions in transcriptional repression, translational regulation and pre-mRNA splicing. For instance, it is reported to regulate alternate splicing of the CTRF gene. In 2006 Neumann et al. found that hyperphosphorylated, ubiquitinated and/or cleaved forms of TDP-43, collectively known as pathological TDP-43, play a major role in the disease mechanisms of ubiquitin-positive, tau- and alpha-synuclein-negative frontotemporal dementia (FTLD-U) and in amyotrophic lateral sclerosis (ALS).

Storage

Storage:
Store at -20°C. Avoid exposure to light.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed mouse brain tissue using CL488-80002 (TDP-43 (for IF/FC) antibody) at dilution of 1:100.



Immunofluorescent analysis of (4% PFA) fixed mouse brain tissue using CL488-80002 (TDP-43 (for IF/FC) antibody) at dilution of 1:100.



1×10^6 HeLa cells were stained with 0.2 μ g CoraLite®488 Anti-Human TDP-43 (for IF/FC) (CL488-80002, Clone:16A22) (red), or 0.2 μ g CoraLite®488 Mouse IgG1 Isotype Control (CL488-66360, Clone: T1F8D3F10) (blue). Cells were fixed and permeabilized with Transcription factor 1X fix/perm buffer.