

For Research Use Only

# CoraLite®594 Anti-Rat CD43 (W3/13)

Catalog Number:CL594-65237



## Basic Information

Catalog Number:

CL594-65237

Size:

100ug , 0.5 mg/ml

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

GeneID (NCBI):

24796

UNIPROT ID:

P13838

Full Name:

sialophorin

Calculated MW:

38 kDa

Purification Method:

N/A

CloneNo.:

W3/13

Recommended Dilutions:

FC: 0.5 ug per 10<sup>6</sup> cells in 100 µl suspension

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

## Applications

Tested Applications:

FC

Species Specificity:

Rat

Positive Controls:

FC : rat splenocytes cells,

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> rat splenocytes cells were surface stained with 0.5 ug CoraLite® 594 Anti-Rat CD43 (CL594-65237, Clone: W3/13), or 0.5 ug Isotype Control. Cells were not fixed.