

For Research Use Only

PE Anti-Human STRO-1 (STRO-1)



Catalog Number: PE-65184 **1 Publications**

Basic Information

Catalog Number:

PE-65184

Size:

100tests , 10 µl/test

Source:

Mouse

Isotype:

Mouse IgM lambda

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

The purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. The conjugate is purified by size-exclusion chromatography.

CloneNo.:

STRO-1

Applications

Tested Applications:

FC

Cited Applications:

IF

Species Specificity:

Human

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yangfan Li	31573683	J Periodontol	IF

Storage

Storage:

Store at 2-8°C. Avoid exposure to light.

Storage Buffer:

TBS with 0.1% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



100 ul human peripheral blood were surface stained with APC Anti-Human CD45 (APC-65109, Clone: HI30) and 10 ul PE Anti-Human STRO-1 (PE-65184, Clone: STRO-1).



100 ul human peripheral blood were surface stained with APC Anti-Human CD45 (APC-65109, Clone: HI30) and PE-Mouse IgM isotype control.