

For Research Use Only.  
Not For Use In Diagnostics.

# Peroxidase-conjugated Affinipure Goat Anti-Horse IgG(H+L)

Catalog Number:SA00001-13



[www.ptglab.com](http://www.ptglab.com)

## Information

Catalog Number:

SA00001-13

Host:

Goat

Applications:

ELISA, WB

Reactivity:

Horse

Physical State:

Liquid

Conjugation:

HRP

## Recommended Dilutions

1:2000-1:10,000 for western blotting with ECL substrates

1:1000-1:20,000 for ELISA and Western blotting with chromogenic substrates

## Stabilizer

3mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

Storage:

Store at -20°C. Avoid repeated freeze / thaw cycles.

Storage Buffer:

0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6

Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data