

For Research Use Only.
Not For Use In Diagnostics.

HRP-conjugated Recombinant Rabbit Anti-Mouse IgG Kappa Light Chain



www.ptglab.com

Catalog Number: **SA00001-19**

Information

Catalog Number:

SA00001-19

Host:

Rabbit

Applications:

ELISA, WB

Reactivity:

Mouse

Physical State:

Liquid

Conjugation:

HRP

Clonality:

Recombinant

Recommended Dilutions

1:2000-1:20000

Preservative

Proclin300

Stabilizer

0.5% BSA

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.03% Proclin300, 50% glycerol, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

Purity

Protein A purification

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data