

For Research Use Only.
Not For Use In Diagnostics.

Alkaline Phosphatase–conjugated Affinipure Goat Anti-Mouse IgG(H+L)



www.ptglab.com

Catalog Number: SA00002-1

Information

Catalog Number: SA00002-1	Reactivity: Mouse
Host: Goat	Physical State: Liquid
Applications: ELISA, WB	Conjugation: Alkaline phosphatase

Recommended Dilutions

1:1000 – 1:3000 for ELISA and Western blotting

Preservative

0.01% Sodium Azide

Stabilizer

5mg/ml BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:
Store at -20°C. Avoid repeated freeze / thaw cycles.
Storage Buffer:
0.01M Tris-HCl, 0.25M NaCl, 50% glycerol, pH8.0
Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data