

For Research Use Only.  
Not For Use In Diagnostics.



[www.ptglab.com](http://www.ptglab.com)

# Goat anti-rabbit IgG (H+L), FITC conjugate

Catalog Number: SA00003-2

## Information

Catalog Number:  
SA00003-2  
Applications:  
Flow Cyt, ICC/IF, IHC

Product name:  
Goat anti-rabbit IgG (H+L), FITC conjugate  
Physical State:  
Liquid  
Conjugation:

## Recommended Dilutions

1:50-1:200 for most applications

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

Storage:  
Store at -20°C. Avoid repeated freeze / thaw cycles.  
Storage Buffer:  
0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6  
Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product, please contact:

T: 1 (888) 4PTCLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.