

For Research Use Only.  
Not For Use In Diagnostics.

# Texas Red-conjugated Affinipure Donkey Anti-Goat IgG(H+L)

Catalog Number:SA00005-3



[www.ptglab.com](http://www.ptglab.com)

## Information

Catalog Number:

SA00005-3

Host:

Donkey

Applications:

IF, FC

Reactivity:

Goat

Physical State:

Liquid

Conjugation:

Texas Red

## Recommended Dilutions

1:20-1:100 for most applications

## Preservative

0.02% Sodium Azide

## Fluorophore

Texas Red sulfonyl chloride (TRSC)  $\lambda_{max}=596nm, \lambda_{em}=620nm$

## Stabilizer

5 mg/ml BSA

## Absorbance Ratio

A596/A280=0.57

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6

Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data