

For Research Use Only.  
Not For Use In Diagnostics.

# R-PE-conjugated Goat Anti-Rabbit IgG(H+L)

Catalog Number:SA00008-2



[www.ptglab.com](http://www.ptglab.com)

---

## Information

Catalog Number:

SA00008-2

Host:

Goat

Applications:

IF, Flow cyt

Reactivity:

Rabbit

Physical State:

Liquid

Conjugation:

R-PE

---

## Recommended Dilutions

1:20-1:100 for most applications

---

## Preservative

0.02% Sodium Azide

---

## Fluorophore

Phycoerythrin Amax=490nm,545nm and 566nm;Emax=580nm

---

## Stabilizer

5mg/ml BSA

---

## Safety Notes

This product is for research use only,not for diagnostic or therapeutic use.

---

## Storage

Storage:

Store at -20°C. Avoid repeated freeze / thaw cycles.

Storage Buffer:

0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6

Aliquoting is unnecessary for -20°C storage

---

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

---

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

---

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data