

For Research Use Only.  
Not For Use In Diagnostics.

# HRP-conjugated Affinipure Goat Anti-Mouse IgA (a)

Catalog Number:SA00012-7



[www.ptglab.com](http://www.ptglab.com)

---

## Information

**Catalog Number:**

SA00012-7

**Host:**

Goat

**Applications:**

ELISA, WB, Dot blot

**Reactivity:**

Mouse

**Physical State:**

Liquid

**Conjugation:**

HRP

**Clonality:**

Polyclonal

---

## Recommended Dilutions

1:250-1:1500 for ELISA

1:1500-1:5,000 for Western blotting with ECL substrate

---

## Stabilizer

3 mg/ml BSA

---

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

---

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 3 mg/ml BSA, pH 7.6.

Aliquoting is unnecessary for -20°C storage

---

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

---

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

---

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data