

For Research Use Only.  
Not For Use In Diagnostics.

# CoraLite594 – conjugated Donkey Anti-Mouse IgG(H+L)

Catalog Number:SA00013-7



[www.ptglab.com](http://www.ptglab.com)

## Information

**Catalog Number:**

SA00013-7

**Host:**

Donkey

**Applications:**

IF, FC

**Reactivity:**

Mouse (Minimal Cross Reactivity:  
Bovine, Chicken, Goat, Guinea Pig,  
Syrian Hamster, Horse, Human,  
Rabbit, Sheep Serum Proteins)

**Physical State:**

Liquid

**Conjugation:**

CoraLite@594

**Clonality:**

Polyclonal

## Recommended Dilutions

1:100 – 1:500

## Preservative

0.02% Sodium Azide

## Fluorophore

CL594 Fluorescent Dye  
A<sub>max</sub>=593 nm, E<sub>max</sub>=614 nm

## Stabilizer

5 mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data