

For Research Use Only.
Not For Use In Diagnostics.

CoraLite Plus 750-conjugated Goat anti-rabbit IgG(H+L)

Catalog Number:SA00015-1



www.ptglab.com

Information

Catalog Number:
SA00015-1
Host:
Goat
Applications:
WB

Reactivity:
Rabbit
Physical State:
Liquid
Conjugation:
CoraLite® Plus 750

Clonality:
Recombinant

Recommended Dilutions

1:5000 – 1:10000

Preservative

0.1% Proclin-300

Fluorophore

CoraLite Plus 750, Amax=755±4 nm; Emax=780±4 nm

Stabilizer

10 mg/mL BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 10 mg/mL BSA, 0.1% Proclin-300, pH 7.4.
Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from culture media supernatant by immunoaffinity chromatography using Protein G beads.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data