

Nur für Forschungszwecke

phospho-Marcks Polyklonaler Antikörper

Katalog-Nr.:10019-1-AP



Allgemeine Informationen

Katalog-Nr.:
10019-1-AP

Größe:
100ul , Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

GenBank-Zugangsnummer:
BC046601

GeneID (NCBI):
17118

Vollständiger Name:
myristoylated alanine rich protein
kinase C substrate

Berechnete Masse:
32 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Maus

Hintergrundinformationen

MARCKS is a member of MARCKS family of PKC substrate and it is also a widely accepted indicator of PKC activation. It shows a ubiquitous tissue distribution. The calculated molecular weight of MARCKS is 39 kDa and its apparent molecular weights were found at 68 kDa to 90 kDa in different species. This protein is a known to regulated cytoskeleton remodeling and membrane recycling. The anti-MARCKS antibody is raised against purified MARCKS protein from mouse brain using a method similar to that previously reported for MacMARCKS

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

