

Nur für Forschungszwecke

TNFAIP1 Polyklonaler Antikörper

Katalog-Nr.: 10096-1-AP

Vorgestelltes Produkt

2 Publikationen



Allgemeine Informationen

Katalog-Nr.: 10096-1-AP	GenBank-Zugangsnummer: BC001643	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 133 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 7126	
Wirt: Kaninchen	Vollständiger Name: tumor necrosis factor, alpha-induced protein 1 (endothelial)	
Isotyp: IgG	Berechnete Masse: 36 kDa	
Immunogen Katalognummer: AG0132		

Anwendungen

Geprüfte Anwendungen:
ELISA

In Publikationen genannte Anwendungen:
IHC, WB

Getestete Reaktivität:
Human, Maus, Ratte

Zitierte Arten:
Human, Maus

Hintergrundinformationen

TNFAIP1 was identified as a human gene activated transcriptionally by tumour necrosis factor alpha in endothelial cells. It is encoded by a 3.5-kilobase transcript and its expression is found to be developmentally regulated in a tissue-specific manner. The open reading frame predicts a 316-amino acid sequence rich in charged residues, particularly at the carboxyl terminus.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Manmeet K Mamik	27726055	J Neuroimmune Pharmacol	IHC
Juliette Sailland	25288732	Proc Natl Acad Sci U S A	WB

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

