

Nur für Forschungszwecke

DAP3 Polyklonaler Antikörper

Katalog-Nr.: 10276-1-AP



Allgemeine Informationen

Katalog-Nr.:
10276-1-AP

Größe:
150ul, Konzentration: 450 µg/ml von
Nanodrop und 260 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG0274

GenBank-Zugangsnummer:
BF220193

GeneID (NCBI):
7818

Vollständiger Name:
death associated protein 3

Berechnete Masse:
46 kDa

Beobachtete Masse:
40-45 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human, Maus

Positivkontrollen:

WB: HeLa-Zellen, humanes Hirngewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



HeLa cells were subjected to SDS PAGE followed by western blot with 10276-1-AP (DAP3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.