

MMP15 Polyklonaler Antikörper

Katalog-Nr.: 10372-1-AP

Allgemeine Informationen

Katalog-Nr.:

10372-1-AP

Größe:

150ul, Konzentration: 247 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG0548

GenBank-Zugangsnummer:

BC006412

GeneID (NCBI):

4324

Vollständiger Name:

matrix metalloproteinase 15
(membrane-inserted)

Berechnete Masse:

72 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

MMP15 (Matrix metalloproteinase-15) is a endopeptidase that degrades various components of the extracellular matrix. MMP15-specific antibodies detect a 72-kD protein in lung cell surface membranes, suggesting that MMP15 is membrane-bound. They demonstrated by Northern blotting that MMP15 is widely expressed as a 3.6-kb transcript, with highest levels in liver, placenta, testis, colon, and intestine (PMID:7649159).

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

