

Nur für Forschungszwecke

Synaptogyrin 1 Polyklonaler Antikörper

Katalog-Nr.: 10432-1-AP



Allgemeine Informationen

Katalog-Nr.:
10432-1-AP

Größe:
150ul, Konzentration: 247 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG0685

GenBank-Zugangsnummer:
BC000731

GeneID (NCBI):
9145

Vollständiger Name:
synaptogyrin 1

Berechnete Masse:
25 kDa

Beobachtete Masse:
28-30 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human, Ratte

Positivkontrollen:

WB: Rattenhirngewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



rat brain tissue were subjected to SDS PAGE followed by western blot with 10432-1-AP (Synaptogyrin 1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.