

Nur für Forschungszwecke

TSPAN3 Polyklonaler Antikörper

Katalog-Nr.: 10447-1-AP

Vorgestelltes Produkt



Allgemeine Informationen

Katalog-Nr.:
10447-1-AP

Größe:
150ul, Konzentration: 247 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG0391

GenBank-Zugangsnummer:
BC000704

GeneID (NCBI):
10099

Vollständiger Name:
tetraspanin 3

Berechnete Masse:
252 aa, 28 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

Tetraspanin-3 (TSPAN3, also known as TM4SF8) is a member of the tetraspanin (TM4SF) family. Tetraspanins, so named because they contain four transmembrane domains, are frequently expressed at the cell surface in association with each other and with other molecules, and they function to regulate cell adhesion, migration, proliferation, and differentiation. TSPAN3 regulates the proliferation and migration of oligodendrocytes, a process essential for normal myelination and repair.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

