

Nur für Forschungszwecke

# PPP1R2 Polyklonaler Antikörper

Katalog-Nr.:10516-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
10516-1-AP

Größe:  
150ul, Konzentration: 200 µg/ml  
durch die Bradford-Methode mit BSA  
als Standard;

Wirt:  
Kaninchen

Isotyp:  
IgG

Immunogen Katalognummer:  
AG0805

GenBank-Zugangsnummer:  
BC007655

GeneID (NCBI):  
5504

Vollständiger Name:  
protein phosphatase 1, regulatory  
(inhibitor) subunit 2

Berechnete Masse:  
23 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

## Anwendungen

Geprüfte Anwendungen:  
ELISA

Getestete Reaktivität:  
Human, Maus, Ratte

## Hintergrundinformationen

PPP1R2 (PHOSPHATASE INHIBITOR 2) forms heterodimer with phosphatase type-1 (PP1). Activity of the complex can regulated by a series of phosphorylations/ phosphorylations events on PPP1R2, functioning in some important process in vivo. The 50kd band recognized by 10516-1-AP seems the complex, but we are not sure about that.

## Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

