

Nur für Forschungszwecke

PLCD4 Polyklonaler Antikörper

Katalog-Nr.:10589-2-AP



Allgemeine Informationen

Katalog-Nr.:
10589-2-AP

Größe:
150ul, Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG0937

GenBank-Zugangsnummer:
BC006355

GeneID (NCBI):
84812

Vollständiger Name:
phospholipase C, delta 4

Berechnete Masse:
85-90 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
IHC, ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunohistochemical analysis of paraffin-embedded human pancreas cancer using 10589-2-AP (PLCD4 antibody) at dilution of 1:50 (under 10x lens).