

Nur für Forschungszwecke

PTGES Polyklonaler Antikörper

Katalog-Nr.:10596-1-AP



Allgemeine Informationen

Katalog-Nr.:
10596-1-AP

Größe:
150ul , Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG0847

GenBank-Zugangsnummer:
BC008280

GeneID (NCBI):
9536

Vollständiger Name:
prostaglandin E synthase

Berechnete Masse:
17 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

Human prostaglandin E synthase (PTGES) is a member of a recently recognized protein superfamily consisting of membrane associated proteins involved in eicosanoid and glutathione metabolism (the MAPEG family) (PMID:10377395). It is also named as MGST1L1, MPGES1, PGES, PIG12 and belongs to the MAPEG family. This protein is expressed in membranes as an approximately 16-kD protein (PMID:10091672).

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

