

Nur für Forschungszwecke

GEMIN7 Polyklonaler Antikörper

Katalog-Nr.:10625-1-AP



Allgemeine Informationen

Katalog-Nr.: 10625-1-AP	GenBank-Zugangsnummer: BC007793	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 750 µg/ml von Nanodrop;	GeneID (NCBI): 79760	
Wirt: Kaninchen	Vollständiger Name: gem (nuclear organelle) associated protein 7	
Isotyp: IgG	Berechnete Masse: 15 kDa	
Immunogen Katalognummer: AG0976		

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 10625-1-AP (GEMIN7 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.