

Nur für Forschungszwecke

TOR3A Polyklonaler Antikörper

Katalog-Nr.: 10647-1-AP

Vorgestelltes Produkt



Allgemeine Informationen

Katalog-Nr.:

10647-1-AP

Größe:

150ul, Konzentration: 500 µg/ml von Nanodrop und 333 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG1012

GenBank-Zugangsnummer:

BC007571

GeneID (NCBI):

64222

Vollständiger Name:

torsin family 3, member A

Berechnete Masse:

46 kDa

Beobachtete Masse:

39-42 kDa, 46-52 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:300-1:1000

IF 1:50-1:500

Anwendungen

Geprüfte Anwendungen:

IF, WB, ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB : HEK-293-Zellen, HepG2-Zellen, PC-3-Zellen, Raji-Zellen

IF : HEK-293-Zellen,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten




HEK-293 cells were subjected to SDS PAGE followed by western blot with 10647-1-AP (TOR3A antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 10647-1-AP (TOR3A antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



WB result of TOR3A antibody (10647-1-AP; 1:2000; incubated at room temperature for 1.5 hours) with sh-Control and sh-TOR3A transfected HEK-293 cells.



Immunofluorescent analysis of (4% PFA) fixed HEK-293 cells using 10647-1-AP (TOR3A antibody) at dilution of 1:50 and Alexa Fluor 488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).